## S.E.R. FAQ NotTaR of Television Sets: CRT purity adjustment

Download

Feedback

Home Copyright © 1994-2007, Samuel M. Goldwasser, All Rights Reserved. Reproduction of this document in whole or in Search part is permitted if both of the following conditions are satisfied: 1. This notice is included in its entirety at the beginning. 2. There is no charge except to cover the costs of copying. I may be contacted via the Sci. Electronics. Repair FAQ (www.repairfaq.org) Email Links Page.

<< CRT purity and convergenc.. | Index | CRT convergence adjustmen.. >>

## **CRT** purity adjustment

Purity on modern CRTs is usually set by a combination of a set of ring magnets just behind the deflection voke on the neck of the CRT and the position of the yoke fore-aft. As always, mark the starting position of all the rings and make sure you are adjusting the correct set if rings!

Use the following purity adjustment procedure as a general guide only. Depending on the particular model TV, the following purity adjustment procedure may substitute green for red depending on the arrangement of the guns in the CRT. This description is based on the Sams' Photofact for the RCA CTC111C chassis which uses a slot-mask CRT. The procedures for dot-mask and Trinitron (aperture grille) CRTs will vary slightly. See you service manual!

Obtain a white raster (sometimes there is a test point that can be grounded to force this). Then, turn down the bias controls for blue and green so that you have a pure red raster. Let the set warm up for a minimum of 15 minutes.

Loosen the deflection voke clamp and move the voke as far back as it will go,

Adjust the purity magnets to center the red vertical raster on the screen.

Move the yoke forward until you have the best overall red purity.

Now, move the yoke forward until you have the best overall red purity. Tighten the clamp securely and reinstall the rubber wedges (if you set has these) to stabilize the voke position. Reset the video adjustments you touched to get a red raster.

5/5/2013 6:14 PM 1 of 1